



extremeH2O[®]

Large (sports) Lens - 14.8

Large Diameter Contact Lenses Have Big Benefits

- Sports Lens For Active Lifestyles

The larger diameter is more stable on the eye and less likely to decenter or fall out during athletic activities.

- Larger Corneas Need More Limbal Clearance

The larger size provides full coverage and sits under the lid. Any lid interaction which could cause discomfort from the lens edge is gone.

Full coverage reduces the feeling of dryness making it a perfect solution for:

- Computer users
- Allergy sufferers
- Middle-aged wearers
- People working in climate controlled offices
- People living in low humidity areas
- Women on birth control and other medications



Why Diameter Matters

Disposable contact lenses generally come in one standard diameter. However, a recent independent study by Pacific University determined that 27% of contact lens fittings fall outside these standard parameters and would benefit from a small or large diameter lens. While a “one size fits all” contact lens is convenient, it overlooks a large segment of your patient base - and might contribute to the steady dropout rate despite strides in contact lens design and breathability.

Ultra-Hydrating Material - Gold Standard for Contact Lens Comfort

Our proven GMA/hydrogel copolymer has advanced hydration properties, enabling it to retain 99% of its moisture - even under extreme conditions. Lens comfort, movement, and fit remain stable all day long.

- Non-ionic material resists deposits keeping lenses clear and free of unhealthy protein deposits
- Naturally hydrophilic - coatings are not required for superior comfort
- Proper DK value ensures good corneal health
- Ultra-stable copolymer provides exceptional on-eye stability throughout the entire wearing period

Extreme H2O 14.8 - Available in Daily and Weekly Replacement

- **Daily:** Premium disposable contact lens
- **Weekly:** Affordable frequent replacement - the next best thing to daily replacement

Product Specifications				
Material	Hioxifilcon A with 59% water			
Polymer Type	GMA/HEMA copolymer, group 2, high water, non-ionic			
Handling Tint	Light blue			
O ₂ Permeability	28 Dk (Fatt units @ 35°C)			
Extreme H2O Daily				
Packaging	Blister pack: 30 per box & 90 per box			
Wear Indication	Daily replacement suggested			
Center Thickness	0.070 mm @ -3.00 mm			
Diameter	13.6 mm	14.2 mm	14.8 mm	
Base Curve	8.2 / Median	8.5 / Median	8.7 / Median	
Powers <i>0.50D steps after -6.00</i>	+0.50 to +6.00 -0.50 to -10.00	+0.50 to +6.00 -0.50 to -10.00	+0.50 to +6.00 -0.50 to -10.00	
Extreme H2O Weekly				
Packaging	Blister pack: 12 per box			
Wear Indication	Daily wear: weekly replacement suggested			
Center Thickness	0.110 mm @ -3.00 mm			
Diameter	13.6 mm	14.2 mm		14.8 mm
Base Curve	8.2 / Median	8.5 / Median	8.3 / Steep	8.7 / Median
Powers <i>0.50D steps after -6.00</i>	+0.50 to +6.00 -0.50 to -10.00	+0.25 to +6.00 -0.25 to -10.00	-0.25 to -10.00	+0.50 to +6.00 -0.25 to -10.00

